



AdVnt Biotechnologies, LLC

17200 W Bell Road, Surprise, AZ 85374

Web Catalog

Plague (*Y. pestis*) Threat Detection Kit 10/bx



Product URL <https://advntbiotechnologies.com/badd-plague-y-pestis-threat-detection-kit-10bx>

Short Description: AdVnt's BADD Plague (*Y. pestis*) Biowarfare Detection Test Kit- 10 Kits/Bx

Full Description: The BADD Plague (*Y. pestis*) Biowarfare Screening Assay is a Department of Homeland Security "Approved Product for Homeland Security" under the Support Anti-terrorism by Fostering Effective Technologies Act of 2002 (the SAFETY Act). This single-agent detection kit is designed to detect *Yersinia pestis* with a limit of detection (LOD) of 100,000 cfu/ml. Each kit contains everything needed to perform the test. No additional collection/dilution kits or readers are required. This item has a 2-year shelf life from date of manufacture and does not require refrigeration. Works in temperatures up to 120° Fahrenheit.

The BADD Plague detection kit includes the following:

- 1- pouched test
- 1- Instruction sheet
- 1- Buffer/Bottle/Dropper tube
- 1- Chain of Custody Label

Price: 300.00 USD

SKU: PLA-KIT-10

