

AdVnt Biotechnologies, LLC

7898 E Acoma Drive, Suite 106 Scottsdale, AZ 85260

Web Catalog

Pro Strips 5 Agent BioWarfare Threat Detection Kit 10/Bx



Product URL <https://advntbiotechnologies.com/pro-strip-5-agent-biowarfare-threat-detection-kit>

Short Description: AdVnt's Pro Strips Rapid Screening System Multi-agent BioWarfare Threat Detection Kit 10/Bx

Full Description: The "Pro Strips-Rapid Screening System" Multi-Agent Biowarfare Screening Assay is a DHS "Approved Product for Homeland Security" under the Support Anti-terrorism by Fostering Effective Technologies Act of 2002 (the SAFETY Act) Pro Strips is the first multi-agent Biowarfare Threat Detection screening Kit that tests for five (5) biological targets (with a single 0.6 mL sample) Anthrax, Ricin Toxin, Botulinum Toxin, Plague and SEB (Staphylococcal enterotoxin B). ThreatLegalLimit of Detection Anthrax Bacillus anthracis; 1.5×10^4 - 8.3×10^8 /mL (strain dependent) Ricin Toxin Ricinus communis 10 ng/mL Botulinum Toxin Clostridium botulinum 33 ng/mL (BoNT A) 500 ng/mL (BoNT B) Plague Yersinia pestis 1.0×10^5 CFU/mL SEB Staphylococcus aureus 10 ng/mL (Bacillus anthracis strains: Vollum, Ames, Sterne & New Hampshire) Each kit contains everything needed to perform the test. No additional Collection/Dilution kits or readers are required. This product has a 2 year (24 months) shelf life from date of manufacture and does not require special storage. Made in the U.S.A./ Pro Strip Rapid Screening System's Commodity Jurisdiction Determination, by Department of State, is not subject to ITAR jurisdiction. The Department of Commerce advises the item is classified as EAR99. Each Pro Strips- Rapid Screening System Kit includes the following:

- 1- Multi-agent test device
- 1- Instruction sheet
- 1- Buffer/swab/dropper tube
- 1- Chain of Custody label

Price: 750.00 USD

SKU: PS-5T

