

AdVnt Biotechnologies, LLC

7898 E Acoma Drive, Suite 106 Scottsdale, AZ 85260

Web Catalog

Pro Strips 5 Agent BioWarfare Threat Detection Kit 10/Bx



Product URL <https://advntbiotechnologies.com/pro-strip-5-agent-biowarfare-threat-detection-kit>

Short Description: AdVnt's Pro Strips Rapid Screening System Multi-agent BioWarfare Threat Detection Kit 10/Bx

Full Description: The "Pro Strips-Rapid Screening System" Multi-Agent Biowarfare Screening Assay is a DHS "Approved Product for Homeland Security" under the Support Anti-terrorism by Fostering Effective Technologies Act of 2002 (the SAFETY Act) Pro Strips is the first multi-agent Biowarfare Threat Detection screening Kit that tests for five (5) biological targets (with a single 0.6 mL sample) Anthrax, Ricin Toxin, Botulinum Toxin, Plague and SEB (Staphylococcal enterotoxin B). ThreatLegalLimit of Detection Anthrax Bacillus anthracis; 1.5×10^4 - 8.3×10^8 /mL (strain dependent) Ricin Toxin Ricinus communis 10 ng/mL Botulinum Toxin Clostridium botulinum 33 ng/mL (BoNT A) 500 ng/mL (BoNT B) Plague Yersinia pestis 1.0×10^5 CFU/mL SEB Staphylococcus aureus 10 ng/mL (Bacillus anthracis strains: Vollum, Ames, Sterne & New Hampshire) Each kit contains everything needed to perform the test. No additional Collection/Dilution kits or readers are required. This product has a 2 year (24 months) shelf life from date of manufacture and does not require special storage. Made in the U.S.A./ ITAR classification is EAR99. Each Pro Strips-Rapid Screening System Kit includes the following:

- 1- Multi-agent test device
- 1- Instruction sheet
- 1- Buffer/swab/dropper tube
- 1- Chain of Custody label

Price: 735.00 USD

SKU: PS-5T

