

AdVnt Biotechnologies, LLC

7898 E Acoma Drive, Suite 106 Scottsdale, AZ 85260

Web Catalog

Pro Strips Advanced Training Kit- ATTAK PAK 10/bx



Product URL <https://advntbiotechnologies.com/pro-strip-attak-pak-professional-training-kit-11bx>

Short Description: Pro Strips ATTAK (Advanced Tactical Threat Assessment Kit) Training System- 10/bx

Full Description: The Pro Strips ATTAK Training system creates the most realistic training experience when preparing to deploy the Pro Strips 5 Agent Biowarfare Threat Detection kits to your biological threat response teams. The ATTAK training kits precisely portrays the use of the Pro Strips Multi-Agent BioWarfare Detection kits. Each test device will provide either positive or negative results for the following BW agents when exposed to the included training simulants (Positive powder, Positive liquid or Negative powder):

- BA- Anthrax
- RC- Ricin toxin
- CB- Botulinum toxin
- YP- Plague(Y. pestis)
- SE- SEB (Staphylococcal enterotoxin B)

Each box includes the following:

- 9 - ATTAK training kits
- 1 - PS-5S Pro Strips Simulation training kit (Shows Positive, Weak positive, Neg & Invalid test results)
- 1 - Positive ThreatDust white simulation powder
- 1 - Positive ThreatJuice simulation liquid
- 1 - Negative ThreatDust white simulation powder
- 1- Trainer instruction guide

Please indicate the pre-programmed test you require by selecting the corresponding product code when ordering PS-A5-SP for random mix of bio-threat outcomes

- PS-A5-BA for Bacillus anthracis
- PS-A5-RC for Ricin Toxin
- PS-A5-CB for Botulinum Toxin
- PS-A5-YP for Y.pestis (Plague)
- PS-A5-SE for SEB (Staphylococcal enterotoxin B)

Price: 580.00 USD

SKU: PS-A5

