



# AdVnt Biotechnologies, LLC

17200 W Bell Road, Surprise, AZ 85374

## Web Catalog

### **Anthrax Threat Detection Kit 10/bx**



Product URL <https://advntbiotechnologies.com/badd-anthrax-threat-detection-kit-10bx>

Short Description: AdVnt's BADD Anthrax Biowarfare Detection Test Kit- 10 Kits/Bx

Full Description: (Reactive to BA strains: Vollum, Ames, Sterne, New Hampshire, Maputo, Turkey and USAMRID)The BADD Anthrax Biowarfare Screening Assay is a Department of Homeland Security "Approved Product for Homeland Security" under the Support Anti-terrorism by Fostering Effective Technologies Act of 2002 (the SAFETY Act)This single-agent detection kit is designed to detect Bacillus anthracis with limit of detection (LOD) that ranges as  $1.55 \times 10^4$ - $1.6 \times 10^8$  cfu/ml (strain dependent)Each kit contains everything needed to perform the test.No additional collection/dilution kits or readers are required.This item has a 2 year shelf life from date of manufacture and does not require refrigeration.Works in temperatures up to 120° Fahrenheit

The BADD Anthrax detection kit includes the following:

- 1- pouched test
- 1- Instruction sheet
- 1- Buffer/Bottle/Dropper tube
- 1- Chain of Custody Label

Price: 300.00 USD

SKU: ANTH-KIT-10

