



AdVnt Biotechnologies, LLC

17200 W Bell Road, Surprise, AZ 85374

Web Catalog

Botulinum Toxin Threat Detection Kit 10/bx



Product URL <https://advntbiotechnologies.com/badd-botulinum-toxin-threat-detection-kit-10bx>

Short Description: AdVnt's BADD Botulinum Toxin Biowarfare Detection Test Kit- 10 Kits/Bx

Full Description: The BADD Botulinum Toxin Biowarfare Screening Assay is a Department of Homeland Security "Approved Product for Homeland Security" under the Support Anti-terrorism by Fostering Effective Technologies Act of 2002 (the SAFETY Act) This single-agent detection kit is designed to detect Clostridium botulinum with limit of detection (LOD) 33 ng/mL (BoNT A), 500 ng/mL (BoNT B) Each kit contains everything needed to perform the test. No additional collection/dilution kits or readers are required. This item has a 2 year shelf life from date of manufacture and does not require refrigeration. Works in temperatures up to 120° Fahrenheit

The BADD Botulinum Toxin Threat Detection Kit includes the following:

- 1- pouched test
- 1- Instruction sheet
- 1- Buffer/Bottle/Dropper tube
- 1- Chain of Custody Label

Price: 300.00 USD

SKU: BOT-KIT-10

