



AdVnt Biotechnologies, LLC

17200 W Bell Road, Surprise, AZ 85374

Web Catalog

Pro Strips Simulation/Trainer Kit 5



Product URL <https://advntbiotechnologies.com/pro-strips-simulation-training-device>

Short Description: Pro Strips Simulation/Trainer Kit 5- 1/kit

Full Description: The Pro Strips Simulation Training device is helpful to demonstrate the potential outcomes of the Pro Strips PS-5T kits. Each kit contains everything needed to perform the test.

This item has a 2 year shelf life from date of manufacture and does not require refrigeration.

Each test contains the following:

- 1 - Test Device
- 1 - Buffer/swab/dropper tube
- 1 - Positive Simulation Control Card
- 1 - Instruction sheet

By swabbing the training card simulant the PS-5S device demonstrates the following examples:
A strong positive result (Test line visible in 3 minutes)
A weak positive result (Faint test line visible in 13 min)
A negative result (Clear test line)
and 2 different invalid results

Price: 60.00 USD

SKU: PS-5S

